

Common Stocks, Antibiotics

Stock, conc	Preparation & Storage (Bennett Lab)	Volumes, Tubes	Location
10x T4 ligase buffer	Thawed, precipitates resuspended with vortexing/slight warming.	50 µL aliquots	Rack 1.5, -20°
10 mM dNTPs		20 µL aliquots	Rack 1.5, -20°
<i>Antibiotics</i>			
chl ²⁵ (mg/mL)	Ethanol, -20° liquid; stable for 1 year.	15 mL	Fr 1 door
amp ¹⁰⁰	40% ethanol, -20° liquid; stable enough for 6 mo ⁽⁸⁾ . Definitely in water, -20° solid for 6 mo ^(10,11) .	15 mL	Fr 1 door
kan ⁵⁰	Water, 0° liquid; stable for 1 year.	15 mL	Fr 0 door
spt ⁵⁰	Water, 0° liquid; stable for 1 year.	15 mL	Fr 0 door
gen ¹⁰	Water, -20° solid; stable for 1 year.	~1 mL	Rack 1.5
str ¹⁰⁰	Water, -20° solid; stable for 1 year.	~1 mL	Rack 1.5, -20°
tet ¹⁰	Ethanol, -20° liquid; stable in dark for 1 year?	~1 mL	Rack 1.5, -20°
zeo ¹⁰⁰	?, -20° solid; stable in dark for 1 year.	~1 mL	Rack 1.5, -20°

Antibiotics

Antibiotic	Class, Target	Abbr*	Stock Concentration		Solvents, Temp	Notes	Rich Medium, common <i>E. coli</i>				Common resistance genes
			<i>m</i> /V	Relative (<i>E. coli</i>)			1× c, Multi-Copy; Range	1× c, Single-Copy	Stability		
Chloramphenicol	misc, ribosome	chl	25 mg/mL	1000x	Ethanol, -20°(l)	ChlR gene is toxic to <i>E. coli</i> in minimal media	25 or 34 mg/L ⁽³⁾	25–340 mg/L ⁽⁵⁾	10 mg/L ⁽¹⁾ or 5 mg/L ⁽⁶⁾	<i>All</i> : 1 yr ^(8–10)	Chloramphenicol acetyltransferase Type A1, <i>cat</i> .
Ampicillin-Na Carbenicillin-2 Na	penicillin, peptidoglycan	amp car	100 mg/mL	1000x	40% ethanol, -20°(l) or water, -20°(s)	40% EtOH keeps it liquid at -20°C. Substitutes carbenicillin.	100 mg/L ⁽³⁾	100–200 mg/L	50 mg/L	<i>Stock</i> : 6 mo with <¼ loss in water ^(10, 11) . <i>Agar</i> : selective for 4 mo ⁽⁸⁾ , 44% loss every 4 mo ⁽⁹⁾ .	-lactamase (<i>bla</i>).
Kanamycin A ·SO ₄	aminoglycoside, ribosome	kan	50 mg/mL	1000x	Water, 0°(l); -20°(s)	Insoluble in 40% ethanol. Precipitates in 35% DMSO.	50 mg/L ⁽³⁾	50–500 mg/L ⁽⁵⁾	20 mg/L ⁽¹⁾	<i>All</i> : 1 yr ^(8–10)	Aminoglycoside 3-phosphotransferase Type I or II, <i>aphA1</i> / <i>aph(3')-I</i> / <i>nptI</i> , <i>aph(3')-II</i> / <i>nptII</i> . Type IIIa works in Gram(+).
Spectinomycin ·2HCl·5H ₂ O	aminoglycoside, ribosome	spt	50 mg/mL	1000x	Water, 0°(l); -20°(s)	Precipitates in 35% DMSO. Using 120 mg/L prevents spontaneous genomic resistance evolving ⁽¹³⁾ .	50 mg/L ⁽³⁾ or 120 mg/L ⁽¹³⁾	50–500 mg/L ⁽⁵⁾	20 mg/L ⁽¹⁾	<i>All</i> : 1 yr ^(8–10)	Aminoglycoside 3-adenylyltransferase, <i>aadA</i> .
Streptomycin ·SO ₄	aminoglycoside, ribosome	str	100 mg/mL	1000x	Water, 0°(l); -20°(s)	Genomic resistance via <i>rpsL</i> mutation is common in lab strains.	100 mg/L			<i>All</i> : 1 yr ^(8–10)	Plasmids: Aminoglycoside 3-adenylyltransferase, <i>aadA</i> . Genomic resistance usually <i>rpsL</i> mutant.
Trimethoprim	diaminopyrimidine, FoIA (DHFR)	tmp	15 mg/mL	1000x	DMSO or 1% acetic acid	Insoluble in water or ethanol. Exogenous thymine can lead to evolution of resistance.	15 mg/L?		10 mg/L ⁽²⁾	?	Dihydrofolate reductase type II, <i>dhfrB2</i> .
Apramycin ·SO ₄	aminoglycoside, ribosome	apr	?	1000x	Water	.	?		?	<i>All</i> : 1 yr ^(8–10)	Aminoglycoside N ³ -acetyltransferase, <i>aac(3)-IV</i> / <i>aacC4</i> .

Gentamicin C-SO ₄	aminoglycoside, ribosome	gen	10 mg/mL?	1000×	Water	.	10 mg/L ⁽⁸⁾		10 mg/L ⁽¹⁾	All: 1 yr ⁽⁸⁻¹⁰⁾	Aminoglycoside-N ³ -acetyltransferase, <i>aac(3)-I</i> / <i>aacC1</i> .
Erythromycin	macrolide, ribosome	ery	50 mg/mL	100×	Water, ethanol, DMSO	.	500 mg/L		200 mg/L ⁽¹²⁾	?	23S rRNA(adenine ²⁰⁸⁵) N ⁶ -methyltransferase, <i>ermC</i> .
Tetracycline-HCl	tetracycline, ribosome	tet	10 mg/mL	1000×	Ethanol	Photolabile.	10 mg/L ⁽³⁾	10–15 mg/L	7.5–15 mg/L ^(1,12)	Stock: 4° 1 wk, -20° 2 mo ⁽¹⁰⁾	metal:tetracycline/H ⁺ antiporter, class C, <i>tetA(C)</i> .
Doxycycline	tetracycline, ribosome	dox			Ethanol	Photolabile?					?
Zeocin (Phleomycin D1-Cu ²⁺)	glycopeptide, DNA	zeo	25/100 mg/mL	1000×	Water	Use 25 mg/L with LB-Lennox (5 g/L NaCl), pH 7.5. Or use 100 mg/L for high-salt LB. Light-sensitive. Limit freeze-thaws.	25 mg/L ^(7, low-salt) 100 mg/L	25–50 mg/L ^(7, low-salt)		Agar: 1 mo (Invivogen)	bleomycin/phleomycin binding protein, <i>ble</i> .
Hygromycin B	aminoglycoside, ribosome	hyg	200 mg/mL	1000×	Water	.	50–200 mg/L ⁽³⁾			?	Hygromycin B phosphotransferase, <i>hph</i> .

Note: aminoglycoside solutions are generally found to be stable for more than a year in 4° water or agar^(8–10).

* Follow [American Society of Microbiology AAC abbreviation guidelines](#)

References

- Poteete, Anthony R., Charles Rosadini, and Christine St. Pierre. "Gentamicin and other cassettes for chromosomal gene replacement in *Escherichia coli*." *Biotechniques* 41.3 (2006): 261-264. <https://doi.org/10.2144/000112242>
- Kittleson, Joshua T., Sherine Cheung, and J Christopher Anderson. "Rapid optimization of gene dosage in *E. coli* using DIAL strains." *Journal of biological engineering* 5.1 (2011): 1-7. <https://doi.org/10.1186/1754-1611-5-10>
- Addgene recommendations <https://www.addgene.org/mol-bio-reference/>
- Qiagen: [Growth Of Bacterial Cultures \(qiagen.com\)](#) Note: these stock and working concentrations are uncommon.
- Sara Molinari, unpublished data, p15A with common genes of *cat*, *aphA1*, *aadA*.
- Marionette Sensor Collection: Meyer, Adam J., et al. "*Escherichia coli* 'Marionette' strains with 12 highly optimized small-molecule sensors." *Nature chemical biology* 15.2 (2019): 196-204. <https://doi.org/10.1038/s41589-018-0168-3>
- Zeocin | Thermo Fisher Scientific
- Bennett Lab experience, where an antibiotic "working" means it kills lab *E. coli* cloning strains in broth, or in agar without forming a lawn of the non-resistant base strain while selecting for colonies harboring the introduced resistance gene.
- Ryan, Kenneth J., et al. "Stability of antibiotics and chemotherapeutics in agar plates." *Applied microbiology* 20.3 (1970): 447-451. <https://doi.org/10.1128/am.20.3.447-451.1970>
- Berendsen, B. J. A., I. J. W. Elbers, and A. A. M. Stolker. "Determination of the stability of antibiotics in matrix and reference solutions using a straightforward procedure applying mass spectrometric detection." *Food Additives & Contaminants: Part A* 28.12 (2011): 1657-1666. <https://doi.org/10.1080/19440049.2011.604045> Sci-Hub.
- Okerman, Lieve, Johan Van Hende, and Lieven De Zutter. "Stability of frozen stock solutions of beta-lactam antibiotics, cephalosporins, tetracyclines and quinolones used in antibiotic residue screening and antibiotic susceptibility testing." *Analytica chimica acta* 586.1-2 (2007): 284-288. <https://doi.org/10.1016/j.aca.2006.10.034>
- Ferrières, Lionel, et al. "Silent mischief: bacteriophage Mu insertions contaminate products of *Escherichia coli* random mutagenesis performed using suicidal transposon delivery plasmids mobilized by broad-host-range RP4 conjugative machinery." *Journal of bacteriology* 192.24 (2010): 6418-6427. <https://doi.org/10.1128/JB.00621-10>
- [Barrick Lab wiki: antibiotics](#)

Lab Improvement Weekly: Investigative Report on Ampicillin. By [Shyam Bhakta](#).

This week's investigative report examines everyone's favorite yet hated bacteriolytic antibiotic: ampicillin.

A 1990 paper¹ discusses how the secreted β -lactamase in a saturated preculture of an Amp strain can quickly consume all the ampicillin in a subculture. A stationary culture of ampicillin resistant cells can have such a concentration of β -lactamase that even a 1/200 – 1/1000 dilution will still contain enough β -lactamase to consume all the fresh ampicillin before all the non-resistant cells from the stationary phase culture have been killed. Even without any initial β -lactamase, enough β -lactamase is secreted quickly that toxic plasmids can quickly be lost midway through growth without facing selection. The authors recommend not allowing cultures to reach stationary phase if you need a high proportion of cells to contain your plasmid. And before subculturing, cells are best washed for Amp selection.

Ampicillin's bactericidal activity has a pH optimum is 5.5–6 and is tenfold less at pH 8³.

A paper on stability of antibiotics in solid medium over time² found that 4 wk-old amp plates stored cold have a 30–50% loss in [amp] (causing a 10% change in zone of inhibition size), BUT no change (especially at our high concentration) that would dip below typical Amp^S strain's MIC. Two-month-old plates would seem to be the limit, as after the resulting fourfold reduction in [amp] would approach the MIC. The rest of our common four antibiotics have no significant changes in agar over 2 months.

Activity reduction of a 100 mg/mL (1000×) aqueous, unbuffered solution of ampicillin³ :

4° 20° -20°

1 d: 15%, 28%

2 d: 33%, 45%

7 d: 65%, 81%

6 mo: <10%

A number of studies have indicated that the stability of ampicillin in solution appears to be a function of pH, temperature, and even the identity of the buffer. Although ampicillin in any form is more readily soluble in base, it rapidly loses activity when stored above pH 7.0. Optimal conditions for storage were suggested as 2–8°C, pH 3.8–5, retaining more than 90% activity for a week. Another review noted that the buffer used can also affect stability: at pH 7, Tris is "highly deleterious to the stability... but not so at pH 5;" citrate is fine at pH 7 but not at pH 5; acetate buffer seems best at pH 6.5 .

Carbenicillin is anecdotally more chemically stable in media, especially acidic and hot, and more resistant to -lactamase, but I can only find published evidence that supports the former ; I'm suspicious/doubt the latter claim. Maybe more carb survives 50° agar, but more amp (I've read 200 mg/L) can be used if it matters. While media stability ought to give carbenicillin an edge in selection against amp-sensitive mutants over long culture periods, especially in acidifying media (glucose carbon source), peptide broths like LB alkalify over culture time, and many-day culturing or a requirement of maintenance of 1× concentration in sterile media is rare. I've found 1× amp in LB30/TB liquid is selective for cloning for at least 2 wk in the fridge, longer for regular LB. The remaining amp only needs to exceed the MIC by a fewfold. 100–200 mg /L is 8–16-fold the MIC.

A new paper deletes a small part of the AmpR RBS so it makes less -lactamase. Lowers how much amp is needed and it reduces -lactamase accumulation to the level that subcultures of AmpR strains from amp to fresh amp medium preserves selection without quickly destroying all the amp .

¹ <https://www.ncbi.nlm.nih.gov/pubmed/2199796?dopt=Abstract>

² <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC376956/pdf/applmicro00107-0177.pdf>

³ <https://www.applichem.com/en/shop/product-detail/as/ampicillin-natriumsalz-ibiochemica/>

Analytical Profiles of Drug Substances, K. Florey, ed., (Academic Press, NY) Vol. 2, 1-61 (1973).

Gallelli, Amer. J. Hosp. Pharm., 24, 425-433 (1967). "Stability studies of drugs used in intravenous solutions, part one."

Lynn, B., Brit. J. Intravenous Therapy, 2, 22 (1981). "The stability and administration of intravenous penicillins."

<https://www.sciencedirect.com/science/article/pii/S0003267006021076?via%3Dihub>

<https://sci-hub.tw/https://doi.org/10.1177/106002807000400802>

<https://doi.org/10.1021/acssynbio.1c00393>