

Oligo Assembly

Oligo assembly and [Polymerase Chain Assembly \(PCA\)](#) are both options for building relatively small dsDNA fragments for use in DNA assemblies. In oligo assembly, oligonucleotides are purchased, annealed, and enzymatically phosphorylated (if not purchased phosphorylated), to produce the fragment. Oligo assembly and PCA are useful in the following cases, where the desired sequence is:

1. Sequence is too small to PCR from a source (50 bp), at which point the fragment be difficult to resolve on a gel from incorrect primer multimer products and can thus be troublesome to purify.
2. Sequence is repetitive or too structured toward the ends where PCR primers would bind. Primers would misprime and make an incorrect PCR product. PCR can still be used to amplify a larger, nonrepetitive segment (as in example 2), so as to save on oligos, but the repetitive segment(s) still must be made by oligo assembly or PCA.
3. Sequence is uneconomical or too slow to obtain from a DNA fragment synthesis service (e.g. Twist, Genscript, IDT gBlocks/synthons). Services seem to have a higher price per base for tiny products. You can make it yourself within hours after receiving your oligos tomorrow.

To be formalized eventually. For now, here's a Slack post from [Shyam Bhakta](#) explaining the process. Somehow I manage to complicate even annealing oligos together.

NEB says to phosphorylate oligos by adding 300 pmol 5 termini with 500 nmol ATP, 1 μ L T4 PNK, 1 \times PNK buffer, in a 50 μ L rxn. PNK has full activity in T4 ligase buffer, which also supplies sufficient ATP. We don't keep ATP aliquots around anyways to use PNK buffer. So just use ligase buffer.

I anneal pairs of oligos *first* to ensure their 5 termini are exposed in the duplex for later phosphorylation. However, in the literature, most people individually phosphorylate oligos first, *assuming* the 5 ends will be exposed in whatever 2° structure the oligo adopts at the 37° phosphorylation temp; then they anneal the phospho-oligos. I don't make this needless assumption, which might invite failure.

Anneal oligo pairs in water, or Tris buffer/TE as I think ought to be better due to buffering without Mg^{2+} , at something like 2.5 μM oligo in something like 10 μL . That's 0.25 μL of each 100 μM oligo (=25 pmol each oligo = 50 pmol total termini). Run the thermocycler anneal program. The long 45 min program is thought to be "better" than the short 20 min one, but in practice, the short one has been fine for making even terminator parts with strong secondary structure. Unnecessary high temp exposure to ssDNA ought to be minimized anyhow. The longest protocol, used classically to anneal DNA probes to a DNA mixture, involves placing the tubes in a beaker of boiling water or in a heated 95°C block, and letting the water bath or heat block cool to room-temp.

Then phosphorylate the annealed oligos by spiking in 10x T4 ligase buffer to a final 1x and then the minimal amount of T4 PNK. 50 pmol termini can be phosphorylated by 1/6 μ L T4 PNK by the previous rule. You can round up to 1/4 μ L since that's the minimal pipettable, or you can make a 2x master mix of buffer and enzyme and add it to the annealing reactions (10 μ L MM + 10 μ L annealed oligos). Incubate 37° 30 min. I let it go for 1–2 hr if I can. Do not heat-inactivate the reaction, which may denature the oligo duplexes.


Now the phosphorylated, annealed ds-oligo product should be at a bit less than the initial 2.5 μM (due to the small added volumes, or 1.25 μM if you added a 2x master mix's volume). A 1:25 dilution into water or TE would put this product just under two-fold ($\sim 100 \text{ nM}$) my prescribed concentration for Golden Gates (50 nM). This can conveniently be done by adding 240 μL water to the top of each rxn. The PCR tube cannot hold any more volume. Any more concentrated the previous reactions' oligos, and you'd need another tube to make a dilution. That or just use $<0.5 \mu\text{L}$ when you use it, or end up using $>25 \text{ fmol}$ in the GG with no issue.

What I thus give you is a *one-pot*, sequential annealing-phosphorylation reaction, more efficient by its use of three less "pots" and tips per rxn and half the enzyme than what my great predecessors taught me.

Annealing:

1. Combine 0.5 μ L each of two 100 μ M oligos with 17 μ L buffer* (e.g. 1.7 μ L 10 \times buffer + 15.3 μ L water).
If using a PNK master mix later, only use 9 μ L buffer.
2. Thermocycle to anneal as described in the tables below.
Or place the tube sin a 95°C block heater and turn the device off to cool slowly to room-temp. Or boil water in a beaker and float the reaction tube in it until the water cools down.

Examples



This promoter with two identical operators at the ends would be impossible to PCR from a template, as each primer has two binding sites in the template and would bind one another and possibly themselves because the operators are roughly palindromic. PCA is also failure-prone because of the same: [Primerize](#) warns of mispriming and suggests two stages of assembly, with separate synthesis of two halves of the promoters which are then combined in a final stage. This may work, but the two 35 bp operator repeats might still cause unexpected mispriming at the final stage due to the large sequence identity in the two fragments.

Oligo assembly is most reliable in this case due to the separate and simple, single-instance annealing of each half of the promoter, completely specified within the oligos – no reliance on DNA polymerase to produce sequence a certain way based on sequence-dependent complexity and specificity of DNA annealing for priming.

The two green oligos and two red oligos are individually annealed in buffer, phosphorylated with PNK, diluted and used directly in the assembly (Golden Gate or plain ligation) as two fragments with exposed ssDNA cohesive/"sticky" ends, aka "overhangs". The outer overhangs connect the fragments to upstream and downstream DNA (in this specific example, the part acceptor vector).



This example shows four fragments, three of which are annealed oligo products (red, orange, blue pairs), and the last fragment is a PCR product (using green primers). The four fragments plus a vector (not shown) were assembled together in a Bsal Golden Gate reaction.

Overhangs A and E connect to the vector, which has its own BsaI sites generating complementary overhangs.

Overhang B connects the two fragments (red and orange) encoding the two halves of the strong hairpin terminator *T917*, which could not be made as a single-pair oligo product, presumably because it has too long of a perfectly palindromic sequence (also making it impossible to prime in PCA or PCR). The two longer oligos perhaps prefer to each form a monomeric hairpin helix than a dimeric helix. Breaking the sequence into two fragments eliminated this problem, as the two fragments split the halves of the palindrome/hairpin.

Overhangs C and D connect the orange and blue oligo fragments with the PCR product generated by the green primers (amplifying from some *DT54* template), after the PCR product's Bsal sites (boxes in green primers) generate the sticky ends C and D in PCR product during the Bsal Golden Gate assembly reaction. *DT54* thus did not have to be made by oligo assembly or PCA, which would require more oligos.

Phosphorylation:

3. Add in 2 μL 10x ligase buffer, then 0.25–0.5 μL T4 polynucleotide kinase (PNK).
Or to 10 μL anneals, add 10 μL PNK master mix: $n \text{ rxns} \times (7.75 \mu\text{L water}, 2 \mu\text{L ligase buffer}, 0.25 \mu\text{L PNK})$
 - If ligase buffer was used for annealing, ATP can be used alone, which gets degraded at high temp. Otherwise just add more ligase buffer to 1x.
4. Incubate 37°C, 30 min. Do not heat-inactivate.¹
5. Add 230 μL buffer* to get a final 250 μL , such that the annealed oligo product is 200 nM.

* $\frac{1}{2}$ –1x DNA elution buffer EB (10 mM Tris-HCl pH 8.5), Tris-EDTA buffer TE, NEBuffer CutSmart or Thermo Tango restriction digest buffer, NEBuffer 3.1 (high-salt), or T4 ligase buffer. See *Annealing Buffer* section.

¹ T4 PNK heat-inactivation isn't useful, and the high temperature risks denaturing the DNA.

Anneal, quick				
	Step	Temp	Time	Notes
	Lid	100°	–	
	Start	94°	1 sec	
	Slow ramp	37°, 0.1°/sec	1 sec	slowly ramp 94° to 20–37°
	End			Move to 4°C and use promptly.

Anneal, long				
	Step	Temp	Time	Notes
	Lid	100°	–	
	Start	94° –0.2°/cycle 0.1°/sec	5–10 sec	–0.1° or –0.2°/cycle, requiring 25 or 50 min, respectively.
	Go to #1	–	285x	$(94^\circ - 37^\circ) \div (0.2^\circ/\text{cyc}) = 285 \text{ cyc}$ More cycles if ramping to lower temp.
	End			Move to 4°C and use promptly.

Thermocyclers can maximally ramp at 0.1°C/sec, which means ramping from 98° to 37° ($T=61^\circ$) would theoretically take a mere 610 sec, 10.2 min. This works for simple anneals. Anneal products that have strong hairpins or repeats or degenerate bases / libraries anneal more accurately if annealed more slowly. This requires cycling down the temp 0.1–0.2°C/cycle (T_{cyc}), holding each cycle for 5–10 sec, starting from full denaturation at 94°C (T_{hi}) and ending 25–37°C (T_{lo}). Calculate the number of cycles: $(T_{\text{hi}} - T_{\text{lo}}) \div T_{\text{cyc}} = n_{\text{cyc}}$.

Annealing Buffer

Annealing oligos in water is a common practice that works. However, it is probably more robust to use a buffer to protect the oligos and exposed cohesive ends, especially when exposed to high temp during annealing. Traditional high NaCl, EDTA annealing buffer is absolutely not necessary or advantageous for oligo assembly, and the salt and EDTA can be bad for downstream usage of the oligo product. Common DNA elution/storage and restriction/ligation buffers make for fine annealing buffers and diluents. They keep oligos at a safe pH and provide helix-stabilizing ions.

For annealing and dilution of oligo products, $\frac{1}{2}$ x concentration may be a bit better when annealing in 10 mM Tris buffer (EB/TE) so that its higher pH (8–8.5) doesn't affect later phosphorylation in pH 7.5 ligase buffer. Ligase buffer is however fivefold the buffering capacity (50 mM), so it doesn't seem to generally matter to use full 1x Tris buffer for simplicity. This of course doesn't matter if the oligos are pre-phosphorylated.

Annealing oligos in ligase buffer or CutSmart/Tango (or another similar restriction buffer) adds in 10 mM Mg^{2+} , which substantially increases the T_m of oligos. NEBuffer 3.1, the high-salt restriction digest buffer, has an additional 100 mM NaCl that also increases the T_m of oligos. Salt (monovalent cations) and the Mg (divalent) help anneal shorter oligos whose T_m would otherwise be too close to room temp or reaction temp (37°C). In such a case, dilute the annealed product with the same 1x NaCl/Mg-containing buffer, as some believe diluting short DNAs with water can denature them.

Very Short Oligo Products / Low T_m

Oligo products are generally designed to be annealed at least several °C above 37° in order to be compatible with phosphorylations and Golden Gate reactions / digestions that operate at 37°. Even higher T_m s would be necessary for products to be used in higher-temp reactions like Gibson isothermal assembly. Oligo products designed with low T_m s, (very short) might anneal at a lower temp, minimally above room temp, 25°. One would probably have to be careful even holding such products in hand, lest it denature them. It might make sense to phosphorylate such oligos *before* annealing, order them phosphorylated, or re-anneal after phosphorylation just to be sure.

Annealing oligos in T4 ligase buffer or CutSmart/Tango (or another similar restriction buffer) adds in 10 mM Mg^{2+} to the reaction, which substantially increases the T_m of oligos. This could help anneal short oligos whose T_m s would otherwise be too close to room temp or reaction temp (37°C) in Tris buffer or water. In such a case, dilute the annealed oligo product with the same Mg-containing buffer.

In other, more typical cases with long duplexes, just plain water is fine for dilution, and even for annealing.